

# Lilly TM-Wick Electrodes

The best choice for non-invasive tympanic ECoChG acquisition

## Electrocochleography at its best

- Brings recording location closer to the ECoChG generator site.
- Non-invasive placement minimizes patient discomfort.
- Recording quality is enhanced as compared to other types of non-invasive electrodes.
- Used for ECoChG recordings for over 20 years.
- Gel tip provides excellent electrode contact in close proximity to the cochlea.
- Thin electrode cable fits alongside foam ear tips when placed in the ear canal.
- Dehydrated electrode for extended shelf life.
- Simple hydration by soaking in saline solution.
- Packaging can be used as a hydrating receptacle.
- Touch-proof electrode lead connector.
- BPA and Latex free.
- Made from biocompatible materials, non-sterile, single patient, prescription use.

Made in the U.S.A. at an ISO13485 certified facility.  
FDA clearance number K073384.

## Ordering information:

- M013910 - Wick Electrode (one electrode).
- M013914 - Wick Cable (one cable).
- M019260 - Wick Electrode Pack (2 cables and 10 electrodes).
- M019261 - Wick Electrodes 50 Pack (50 electrodes)
- M019262 - Wick Electrodes 100 Pack (100 electrodes)



Phone: +1 (305) 668-6102 • Toll Free: 1 (800) 447-9783 (USA & Canada)  
Fax: +1 (305) 668-6103 • E-mail: [sales@ihsys.com](mailto:sales@ihsys.com) • Web: [www.ihsys.com](http://www.ihsys.com)

**CE**  
Product Conforms with  
93/42/EEC MDD  
TMW-BRCH-REV-55799