

Video OTOSCOPY

Let Our Image Enhance Your Image



Superior Image Quality
Is The Best Way To Describe
Our Video Otoscope
System

MedRx[®]

www.medrx-usa.com or Call Toll Free: 888-392-1234

Video Companion

Bringing Video Otoscopy To Your Computer



VO with External Light source shown above

The remarkably clear resolution is particularly helpful in detecting excessive cerumen, foreign bodies, external otitis, PE tube status, cholesteatoma, TM perforations and performing cerumen management

- **Deluxe 3mm Otoscope Probe**
- **Compact Color Video Camera**
- **Focus Always Feature**
- **0-150 Watt Halogen Light Source**
- **USB 2.0 Video Capture Cable**
- **NOAH™ & TIMS® Compatible**

The alternative battery operated LED light source is designed for the practitioner that requires portability



MedRxINC.®

www.medrx-usa.com

1200 Starkey Rd. #105,

Largo, FL 33771

Toll Free 888-392-1234

727-584-9600 • Fax: 727-584-9602

CE

OS-S-ADH093-1

MedRx's Authorized Representative in Europe
(Regulatory affairs only) Emergo Europe • Molenstraat 15
2513 BH The Hague • The Netherlands
Tel: (31) (0) 70 345-8570 • Fax: (31) (0) 70 346-7299